Xblitz®

USER'S MANUAL [ALCONTROL MINI]

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User's manual in English language - EN

01. PRODUCT DESCRIPTION

1. Mouthpiece
2. Screen
3. Battery lid (on the back of the device)
4. Air outlet (do not cover)
5. Power on / off

NEVER DRIVE AFTER ALCOHOL CONSUMPTION

02. SCREEN

READY

WARNING

LOW BATTERY
LEVEL

EN

1. READY

Green diode turns on when the device is ready to run a test. Diode will flash if the test result is below warning level.

2. LOW BATTERY LEVEL

Yellow diode starts flashing if the battery voltage is too low. Insert new alkaline batteries 1,5V.

3. WARNING

If the test result shows more than 0.25mg/l / 0.50‰, red warning diode will flash.

4. Measurement unit is set depending on the region.


5. Test result is expressed in mg/l and per mills.

03. USAGE

1. Insert mouthpiece into an adequate inlet.
2. Press ON/OFF button located underneath the screen. Firstly, screen will be turned on and „on” command will be displayed together with welcome, sound signal. Subsequently, device starts counting that lasts 15-30 seconds. This heat-up process prepares the sensor and electronics of the device to run the test.
3. When you hear „beep” sound, and „bLo” command will be displayed start constant blowing into the mouthpiece until you hear „beep” sound again and finish blowing. At this point, if you do not blow air for 20 seconds, device will turn itself off automatically. Turn it back on to run another test.
4. For 3-22 seconds there will be spinning circle icon shown on the screen, and then test result will be displayed together with „beep” sound. If the test result is above 0.25mg/l / 0.50‰ (default settings) red warning diode will be flashing and alert sound will be on. Test result will be visible for 20 seconds (first 15 seconds green or red diode will flash, last 5 seconds diode will be off). If you press „Power” button when the green or red diode is off, device will go back to counting state before another test.
5. Finally, device will display „OFF” command and turns itself off automatically together with adequate sound.
6. To run the test again, repeat steps: 2-4.
7. Interval between tests should equal at least 2 minutes. After finishing test turn the device off by pressing „power” button for at least 2 seconds.

04. BATTERY REPLACEMENT

Remove ribbon from the battery before first use.


Use little power of pressure on the battery lid located on the back of the device, then press (1) and eject it (2) from the device to unlock the lid.

EN


Place alkaline batteries AAA in accordance with below pattern.

Step 1


Insert battery starting from the negative (-) side.


Step 2


Step 3

Step 4

Step 5

NOTE: Make sure that the battery is in a correct positive (+) and negative (-) position, just like shown on a picture.

05. SAFETY ISSUES

- 1.** Wait at least 20 minutes after eating or drinking before you run a test. All kinds of alcohol remaining in mouth or saliva may interfere correct measurement and need ~20 min to dissipate. This also refer to different products made on the basis of alcohol such as mouthwash liquids and others.
- 2.** Avoid making tests in conditions of strong wind or closed room with polluted air.
- 3.** When „BAT” command is displayed it is recommended to replace alkaline batteries 1,5V.
- 4.** Do not blow cigarette smoke, food or liquids into the device, as it may cause damage or destruction of the sensor.
- 5.** Do not use this device as a tool defining whether you or any other user should operate any vehicle or machine.
- 6.** Store device in temperature range: 0 – 40 °C. Device is made to operate within temperature range: 5 – 40 °C. Store and use the device within recommended temperature scopes.
- 7.** Correlation between alcohol content in an exhaled air and alcohol concentration in blood depends on many variables, including environment factors (such as air quality, wind, humidity, temperature etc.) as well as health condition of the person tested.
- 8.** At normal usage device demands renewed calibration once every 6 months (not less than once per 12 months) or after 500 tests, in order to obtain accurate results. Contact your seller or service centre to renew calibration of your device.
- 9.** Do not use test results from this device in front of court. Driver bears full responsibility for driving under influence of alcohol.
- 10.** Do not use this device in combination with ignition blockade in vehicle.
- 11.** Test result may differ depending on many factors. Obtaining given result at once does not mean that the result will be exactly the same few minutes later.

12. Do not let yourself or anyone else drive after alcohol consumption. Always use taxi or other vehicle with sober driver.

NOTE: Never blow directly after alcohol consumption. Such action may cause burning of the sensor and is not covered by warranty.

06. TECHNICAL SPECIFICATION

- Sensor:** Advanced sensor based on fuel cell
- Measurement accuracy:** $\pm 0.005\%$ BAC at 0.050% BAC
- Detection range:** $0.00 \sim 2.00\text{mg/l BrAC}$ / $0.00 \sim 4.00\text{‰}$
- Heat-up time:** $15 \sim 35$ seconds
- Reaction time:** $3 \sim 22$ seconds
- Dimensions:** $111\text{mm} \times 35\text{mm} \times 24\text{mm}$
- Weight:** 51.1g , (73.7g with battery)
- Power supply:** Alkaline batteries 1.5V "AAA" x2
- Working temperature:** $5^\circ \sim 40^\circ$
- Storage temperature:** $0^\circ \sim 40^\circ$
- Calibration:** Each 12 months (we recommend each 6 months) or after 500 tests
- Warranty:** 1 year of limited warranty

07. WARRANTY / CLAIMS

The product is covered by the 60-months warranty. Warranty conditions can be found on the website: <https://xblitz.pl/warranty>

Complaints should be reported via the complaint form located on the website: <https://reklamacje.kgktrade.pl/>

Details, contact and service address may be found at: www.xblitz.pl

Serial number.....

The warranty is only valid for regular servicing and calibration of the device at an authorized service center. No less than once a year. The device from the novelty is properly calibrated. Each breathalyzer offered on the market requires periodic calibration, ie every 12 months (recommended every 6 months), the device should be configured to correctly interpret the amount of alcohol in the exhaled air. The need for recalibration is not a device disadvantage and is not covered by the warranty. The device not regularly calibrated will show incorrect results.

EN

Importer:

KGK Trade Sp. z o.o. sp.k
os. Urocze 12
31-953 Cracow
Poland
Made in South Korea

RoHS

CALIBRATION

<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>
<i>Seller's details:</i>	<i>Calibration date:</i>

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